


# TLA



## High Power 7 Colour Neck and Face Mask

### Short Description

This Photon Beauty Device, is designed for facial beautification and uses natural light wavelengths, which are transmitted by LEDs into the skin. The light activates photoreceptors in the skin cells, producing energy for absorption by skin components, helping to beautify your complexion.

### Description

Using the latest in LED technology, we have developed this LED mask to be both comfortable and functional. We recommended that you cleanse your face and neck, and apply a non-light sensitive mask or moisturiser before use.

Over the last few years, the use of daily/ weekly light therapy treatments and LED masks like this have become common in the beauty industry. Our Highly Powered 7 colour LED face and neck mask makes it easier and more convenient to experience facial light therapy at home. This models' light weight design makes for comfortable and effective home usage, and could save you \$1000s on beauty salon treatments.

When using this product, we recommend testing each colour for a short period of time, to make sure your skin is not sensitive to any of these wavelengths of light. We also recommend light blocking goggles if you find light is penetrating through to your eyes.

## User Guide

### Follow these steps

1. Plug your mask in to the wall and its remote.
2. Place the mask on in a comfortable position.
3. Turn power on and select colour using your remote
3. Test each wavelength/ colour on your skin for a short time, to make sure no adverse reaction occurs. To change colours, power output, or time settings, simply use the remote
4. For best results: cleanse skin, and apply a mask before using this product, then simply set the light to the colour you want to use. It is beneficial to apply a moisturiser first, and make sure you don't use any cosmetics that cause light sensitivity BEFORE use! Save them for after.
6. Use your remote to turn your product off, and unplug between uses
7. Don't bend the mask in any way. Be very gentle with it and store it in a safe place.

## A Guide to the 7 Colours

**Blue** – Helps to eliminate bacteria

**Red** – May lead to an increasing in collagen production

**Yellow** – alleviates redness from skin

**Green** – Balances the skin by reducing pigmentation

**Orange** – Revitalises and boosts the skin

**Cyan** – Soothes irritated skin

**Purple** – Increases facial rejuvenation

## Time Guide

You can also use your mask up to 4 times a week for 20 minutes.

Use for 21 days minimum for best results. After 21 days, your skin should definitely see the results! Everyone's body reacts differently to different lights, so experiment and test to find the best colours for your body!

## **Suitable for people:**

With an irregular light diet or sleeping patterns, or frequent tech users (office workers, social media experts, etc.)

Not suitable for people with sensitivities to light.